

DV-E Viscometer

Viscosity by Brookfield
the world standard



Flip to "Speed"
Turn the knob
Choose RPM – Done

Displayed info:
Viscosity (cP or mPa•s)
% Torque
Speed/Spindle

Spindle/Speed
Selection

Flip to "Spindle"
Turn the knob
Choose Spindle

Auto Range
Push for determining
highest range capability

Think

MODEL	VISCOSITY RANGE* cP(mPa•s)		SPEEDS	
	Min.	Max.	RPM	Number of Increments
LVDVE	1 [†]	2M	.3-100	18
RV DVE	100 ^{††}	13M	.3-100	18
HADVE	200 ^{††}	26M	.3-100	18
HBDVE	800 ^{††}	104M	.3-100	18

* Using spindles provided. Ranges can be expanded or reduced using accessories.
† 1 cP achieved with UL Adapter accessory. 15 cP on LV with standard spindles.
†† Minimum viscosity is achieved with optional RV/HA/HB-1 spindle.

